

## Jet Massage

Flexible massage station in several different varieties

- · 2 different nozzles, soft and medium
- Combine as you like with 1, 2, 4 or 6 nozzles in a massage station
- · For public and private pools
- Meets the requirements for pool standards
   EN 16582 and EN 13451
- · Available both in Marine and Classic

NEW!
A luxurious
massage
station
for the pool









Marine Jet Massage



Products from Pahlén Pool System, for a trouble-free pool ownership.

For more information: www.pahlen.com



## Upgrade your pool with Pahlén Jet Massage

Jet Massage consists of a standalone inlet, outlet and start button that combines as desired and up to 6 nozzles. All nozzles are made of acid-proof stainless steel AISI 316L and is available in both Marine and Classic designs. The massage station can be controlled with a Pahlén circulation pump (e.g. FloVISE or P2000) and Pahlén control box.

| Item no.   | Туре  | Wood/plastic/<br>steel with liner | Concrete<br>liner | Pipe<br>connection                        |
|--|---|-----------------------------------|-------------------|---|
| 13005571   | Start button Classic 350  |                                   | Х                 |   |
| 13005570   | Massage Soft Classic 350  |                                   | х                 | Ø50mm                                     |
| 13005572   | Massage Medium Classic 350  |                                   | X                 | Ø50mm                                     |
| 13005561   | Start button Classic  | x                                 |                   |   |
| 13005560   | Massage Soft Classic  | x                                 |                   | Ø50mm                                     |
| 13005562   | Massage Medium Classic  | x                                 |                   | Ø50mm                                     |
| 13005551   | Start button Marine 350   |                                   | Х                 |   |
| 13005550   | Massage Soft Marine 350   |                                   | х                 | Ø50mm                                     |
| 13005552   | Massage Medium Marine 350   |                                   | Х                 | Ø50mm                                     |
| 13005541   | Start button Marine   | x                                 |                   |   |
| 13005540   | Massage Soft Marine   | х                                 |                   | Ø50mm                                     |
| 13005542   | Massage Medium Marine   | x                                 |                   | Ø50mm                                     |
| 13005572<br>13005561<br>13005560<br>13005562<br>13005551<br>13005550<br>13005552<br>13005541<br>13005540 | Massage Medium Classic 350 Start button Classic Massage Soft Classic Massage Medium Classic Start button Marine 350 Massage Soft Marine 350 Massage Medium Marine 350 Start button Marine Massage Soft Marine | x<br>x<br>x                       | x<br>x<br>x       | Ø50mn<br>Ø50mn<br>Ø50mn<br>Ø50mn<br>Ø50mn |



Jet Massage in a Classic model

636200 Massage nozzle Soft as option 636201 Massage nozzle Medium as option

## Massage station configurations

Complete Jet Massage consists of 1-6 massage nozzles, a start button, at least 2 outlets, a separate circulation pump, a control box, connections and pipes. The number of massage nozzles governs the size of the pump and the number of outlets.

| Number of massage nozzles | Number of outlets | Pump            |
|---------------------------|-------------------|-----------------|
| 1                         | 2                 | FloVISE 0,55 kW |
| 2                         | 2                 | FloVISE 1,1 kW  |
| 4                         | 3                 | P2000 2,2 kW    |
| 6                         | 4                 | P2000 4,0 kW    |

